

Investigating the relationship between the concentration of MnO_2 catalyst powder and the rate constant in the decomposition of H_2O_2 .

Research Question: How does changing the concentration of solid MnO_2 catalyst powder suspended in aqueous 2% (v/v) H_2O_2 (aq) affect the first-order rate constant for its decomposition at room temperature?

August 18, 2025

Word count: 123

Materials:

- 2 % (v/v) H_2O_2 (aq)
- $\text{MnO}_2(s)$ powder, about 0.05, 0.1, 0.15, 0.2 g

Using Vernier Pressure sensor¹

1 Introduction

Hydrogen peroxide is widely used in disinfection, bleaching, and environmental remediation. Its catalytic decomposition improves its oxidizing ability and is especially relevant in green chemistry due to the reduced waste generated. Transition metal ions such as manganese, iron, cobalt, and lead are commonly used to catalyze this reaction.[?] Manganese, in the form of MnO_2 , is the catalyst of choice in this investigation because it is available in the school laboratory, being stable, inexpensive, and non-toxic in small quantities.

2 Background Information

Most existing studies use MnO_2 catalysts immobilized on solid supports[?] or embedded in porous materials[?]. In contrast, this investigation uses fine MnO_2 powder suspended in solution. This setup introduces practical limitations: the catalyst cannot be recovered after the reaction, and the surface area in contact with the solution is difficult to estimate precisely.

Nonetheless, the powdered form offers simplicity for a school laboratory. Since

¹<https://www.vernier.com/product/gas-pressure-sensor/> (archived), Accessed: 2025-07-04

previous research has shown that the rate constant for MnO_2 -catalyzed H_2O_2 decomposition does not strongly depend on temperature, temperature control will not be a central concern in this investigation.

Existing studies monitor the reaction utilizing methods such as titration on samples taken at specific intervals[?] or colorimetric methods[?], both of which allow direct measurement of the concentration of hydrogen peroxide in the solution. Instead, this experiment opts to evaluate the progress of the reaction by measuring the amount of (oxygen) gas released by the reaction over time. This can be calculated using the ideal gas law given constant temperature and volume and information about the pressure within the system over time. Gas pressure will be measured using a pressure sensor, allowing for precise, continuous data collection suitable for computer processing. This method should require less labor and resources compared to either method mentioned above, since there is no need for additional reagents or equipment that may not be available at my school, nor is there any for manual input while the reaction is taking place.

2.1 Relevant Theory


“The rate limiting reaction is believed to be the initial reaction between hydrogen peroxide and the iron oxides [?, ?, ?]. Additionally, the results reveal that the decomposition of H_2O_2 follows pseudo-first order kinetics:

$$-\frac{d[\text{H}_2\text{O}_2]}{dt} = k_{\text{app}}[\text{H}_2\text{O}_2] \quad (2.1)$$

[...] where k_{app} is the apparent first-order rate constant, and $[\text{H}_2\text{O}_2]$ is the concentration of H_2O_2 in the solution at any time t .”[?]

Let $n(\text{H}_2\text{O}_2)$ and $n(\text{O}_2)$ denote the amount of H_2O_2 present in the solution and the amount of O_2 gas in the flask released by the reaction, respectively. The catalysed decomposition

reaction taking place is described by



If n_0 is the initial amount of H_2O_2 in the solution, we have

$$\begin{aligned} n(\text{O}_2) &= 2(n_0 - n(\text{H}_2\text{O}_2)) \\ \implies n(\text{H}_2\text{O}_2) &= n_0 - \frac{n(\text{O}_2)}{2} \end{aligned} \tag{2.2}$$

Let V_l and V_f denote the solution volume and the flask's total volume, respectively. By the definition of concentration, we have

$$[\text{H}_2\text{O}_2] = n(\text{H}_2\text{O}_2) \cdot V_l^{-1} \tag{2.3}$$

Let P and T denote the pressure and temperature of gas within the flask, respectively. If R is the ideal gas constant, the ideal gas law states

$$n(\text{O}_2) = P(V_f - V_l) \cdot (RT)^{-1} \tag{2.4}$$

Our goal is to compute an approximate value of the constant k_{app} given the data of pressure over time. Equation (2.1) shall be algebraically manipulated with the aid of equations (2.2), (2.3), and (2.4), to express k_{app} in terms of known constants and our function $P(t)$. The following calculations assume that T and V_l are approximately constant over time. ²

² V_f is obviously constant over time within a trial, but not between trials because I seal the flask with a cork.

$$k_{\text{app}}[\text{H}_2\text{O}_2] = -\frac{d}{dt}[\text{H}_2\text{O}_2] \quad (2.1)$$

$$k_{\text{app}} = -[\text{H}_2\text{O}_2]^{-1} \cdot \frac{d}{dt}[\text{H}_2\text{O}_2] \quad (2.3)$$

$$= -(n(\text{H}_2\text{O}_2) \cdot V_l^{-1})^{-1} \cdot \frac{d}{dt}n(\text{H}_2\text{O}_2) \cdot V_l^{-1}$$

$$= -V_l \cdot n(\text{H}_2\text{O}_2)^{-1} \cdot V_l^{-1} \frac{d}{dt}n(\text{H}_2\text{O}_2)$$

$$= -n(\text{H}_2\text{O}_2)^{-1} \cdot \frac{d}{dt}n(\text{H}_2\text{O}_2)$$

$$= -(n_0 - \frac{n(\text{O}_2)}{2})^{-1} \cdot \frac{d}{dt}(n_0 - \frac{n(\text{O}_2)}{2}) \quad (2.2)$$

$$= -(n_0 - \frac{n(\text{O}_2)}{2})^{-1} \cdot (-\frac{1}{2}) \frac{d}{dt}n(\text{O}_2)$$

$$= \frac{1}{2}(n_0 - \frac{n(\text{O}_2)}{2})^{-1} \cdot \frac{d}{dt}n(\text{O}_2)$$

$$= \frac{1}{2}(n_0 - \frac{P(V_f - V_l) \cdot (RT)^{-1}}{2})^{-1} \cdot \frac{d}{dt}P(V_f - V_l) \cdot (RT)^{-1} \quad (2.4)$$

$$= \frac{1}{2}(n_0 - \frac{P(V_f - V_l) \cdot (RT)^{-1}}{2})^{-1} \cdot (V_f - V_l) \cdot (RT)^{-1} \frac{d}{dt}P$$

$$= \frac{1}{2}(n_0 - \frac{P(V_f - V_l) \cdot (RT)^{-1}}{2})^{-1} \cdot \frac{V_f - V_l}{RT} \cdot \frac{dP}{dt}$$

$$= (2n_0 - P(V_f - V_l) \cdot (RT)^{-1})^{-1} \cdot \frac{V_f - V_l}{RT} \cdot \frac{dP}{dt}$$

$$= \frac{V_f - V_l}{RT(2n_0 - P(V_f - V_l) \cdot (RT)^{-1})} \cdot \frac{dP}{dt}$$

$$= \frac{V_f - V_l}{2n_0RT - (V_f - V_l)P} \cdot \frac{dP}{dt}$$

Thanks to the ideal gas law we can substitute the expression $P_0(V_f - V_l) = n_0RT$ where P_0 is the initial pressure of the gas within the flask, to simplify and obtain

$$k_{\text{app}} = \frac{1}{2P_0 - P} \cdot \frac{dP}{dt} \quad (2.5)$$

3 Measurements

$$V_f = 0.345l$$

$$V_l = 0.050l$$

Addressing random variance in measurements: Due to random experimental fluctuations, the raw pressure data P exhibits noise which causes the numerical derivative

$\frac{dP}{dt}$ to occasionally take unphysical negative values. To reduce noise and obtain physically meaningful estimates of the rate of pressure change $\frac{dP}{dt}$, I applied a five-point Savitzky-Golay filter of degree 2 using fixed convolution weights. This method smooths the data while computing derivatives more accurately than basic finite difference formulas. I used the coefficients $\frac{[-2, -1, 0, +1, +2]}{10\Delta t}$, which are derived from a least-squares fit of a quadratic polynomial over a five-point window.