



***Uncertainties and  
errors in  
measurements and  
results***

# Scientific research

*When a measurement is recorded, there is always an experimental error or random uncertainty associated with the value.*

**Uncertainty** is the quantification of the doubt in measurement results.

**Error** can be random or systematic.






## GCSE Working Scientifically "Uncertainty"

# Estimating the uncertainty

- $\pm$  half of the smallest division on the scale
- $\pm 1$  the last significant figure in a digital measurement
- Check data provided by the manufacturer

***Remember to use SI units!***






**51.4 Figure 2** Glassware commonly used in the laboratory. (a) conical or Erlenmeyer flask, its shape makes it easy to mix liquids as the flask can be easily swirled; (b) beaker; (c) measuring or graduated cylinder; (d) volumetric flask; (e) pipette; (f) burette.

Glassware	Volume / $\text{cm}^3$	$\pm$ Uncertainty / $\text{cm}^3$	Uncertainty / %
beaker	50	5	10
measuring cylinder	50.0	0.5	1
burette	50.00	0.05	0.1

## Example 1.

*What's the value and measurement uncertainty in this case?*


**$24.20 \pm 0.05$  ml**



## Example 2.

*What's the value and measurement uncertainty in this case?*

**$39.0 \pm 0.5$  ml**



## Example 3.

*What's the value and measurement uncertainty in this case?*

$5.67 \pm 0.01 \text{ g}$





# How to evaluate the accuracy of your results?

→ **Percentage error:**

$$\text{percentage error} = \left| \frac{\text{literature value} - \text{experimental value}}{\text{literature value}} \right| \cdot 100\%$$

→ **Percentage error tells you how close you were to the correct value.**

# Uncertainty propagation based on the measurement uncertainties

## → Absolute uncertainty:

*For example : initial temperature =  $34.55 \pm 0.05 \text{ }^{\circ}\text{C}$*

## → Relative uncertainty:

$$\text{Relative uncertainty} = \frac{\text{Absolute uncertainty}}{\text{Magnitude of the measurement}}$$

## → Percentage uncertainty:

$$\text{Percentage uncertainty} = \frac{\text{Absolute uncertainty}}{\text{Magnitude of the measurement}} \cdot 100 \%$$

# Propagation of uncertainty

Uncertainty propagation follows your calculations in the data processing:

**If you Add and/or subtract in the data processing:**

- Add the absolute uncertainties

**If you Multiply and/or divide in the data processing:**

- Convert to relative / percentage uncertainty
- Add the relative / percentage uncertainties
- Convert back to absolute uncertainty



## Experimental error based on trials

- “*Uncertainty for average values*”

→ Experimental error: 
$$\frac{\text{max} - \text{min}}{2}$$

→ Experimental error tells about the precision of your results.

# Accuracy and precision

Accuracy = How close a measurement is to the true value

Precision = How close replicate measurements are to each other



✓ Precision  
✗ Accuracy



✗ Precision  
✓ Accuracy



✗ Precision  
✗ Accuracy



✓ Precision  
✓ Accuracy

## To reduce uncertainty

- When designing an experiment, choose equipment that reduces the level of uncertainty
  - e.g. measure 100 ml solution with graduated cylinder, not with beaker
- Do experiment many times
  - taking more readings and calculating the average

## Example 4.

The following data are collected during a titration.

$$\text{Final burette reading} = 16.10 \pm 0.05 \text{ cm}^3$$

$$\text{Initial burette reading} = 1.10 \pm 0.05 \text{ cm}^3$$

Calculate the percentage uncertainty of the titre.

$$\text{percentage uncertainty} = \frac{\text{uncertainty}}{\text{value}} \cdot 100 \%$$

[1

## Example 4.

The following data are collected during a titration.

$$\text{Final burette reading} = 16.10 \pm 0.05 \text{ cm}^3$$

$$\text{Initial burette reading} = 1.10 \pm 0.05 \text{ cm}^3$$

Calculate the percentage uncertainty of the titre.

$$\text{percentage uncertainty} = \frac{\text{uncertainty}}{\text{value}} \cdot 100 \%$$

$$\llcorner \frac{0.1 \text{ cm}^3}{15.0 \text{ cm}^3} \times 100 = \rceil 0.7 \llcorner \% \rceil \checkmark$$



## Example 5.

What is the uncertainty of the concentration of the solution with 2.70 g anhydrous iron(III) chloride, when the uncertainty of the scale is  $\pm 0.01$  g and the uncertainty of the volume is  $\pm 0.001$  dm<sup>3</sup>? (2.70 grams of the solid compound is dissolved to distilled water and diluted up to 0.500 dm<sup>3</sup> in a volumetric flask)

$$c_{\text{anhydrous}}(\text{FeCl}_3) = \frac{n(\text{FeCl}_3)}{V_{\text{solution}}} = \frac{0.016646 \text{ mol}}{0.500 \text{ dm}^3} = 0.0333 \text{ mol dm}^{-3}$$

percentage uncertainty of the mass

$$\frac{0.01 \text{ g}}{2.70 \text{ g}} \times 100\% = 0.37 \%$$

percentage uncertainty of the volume

$$\frac{0.001 \text{ dm}^3}{0.500 \text{ dm}^3} \times 100\% = 0.20 \%$$

percentage uncertainty of the concentration

$$0.37 \% + 0.20 \% = 0.57 \%$$

Absolute uncertainty of the anhydrous solution is then

$$\frac{0.0333 \text{ mol dm}^{-3}}{100 \%} \times 0.57 \% = 0.00018977, \text{ thus}$$

$$c_{\text{anhydrous}}(\text{FeCl}_3) = 0.0333 \text{ mol dm}^{-3} \pm 0.0002 \text{ mol dm}^{-3}$$